


Stable Isotope Gas Analyzer

Continuous measurement of ^{13}C and ^{18}O isotopes in exhaled breath

The new Stable Isotope Gas Analyzer is a cutting-edge upgrade to our standard Promethion™ systems, allowing simultaneous measurement of stable isotope tracers synchronously with the Promethion data stream.

Measure the oxidation of *exogenous* nutrients

Figure A shows the oxidative disposal of a 2.5 mg bolus of ^{13}C -glucose in a control mouse. Possible experimental manipulations include: age, diet, hormonal treatment, experimental drugs, microbiome manipulations, exercise, thermal exposure, surgical procedures, illness/injury, etc. Critical metrics include magnitude and duration of response and AUC for calculating % dose recovery.



FEATURES

Simultaneous measurement of $^{13}\text{CO}_2$, C^{18}O_2 , CO_2 , and H_2O

One analyzer can be multiplexed with up to 8 cages

Fast response time and low power requirement

Wide measurement ranges for $\delta^{13}\text{C}$ (-100‰ to 4000‰) and CO_2 (300 ppm to 25,000 ppm)

High precision for both $\delta^{13}\text{C}$ (0.15‰) and $\delta^{18}\text{O}$ (1.0‰) as well as CO_2 (0.05 ppm) and H_2O (50 ppm)


Over 10,000 different isotopically-labeled tracers commercially available

No consumables or external calibration required

Stable Isotope Gas Analyzer

Measure the oxidation of *endogenous* nutrients

Figure B shows the oxidation of a fatty acid tracer (^{13}C -Palmitic acid) infused into rodent diet for 10 days, selectively enriching the body lipids with ^{13}C . The CO_2 and $\delta^{13}\text{C}$ measured during rest, treadmill exercise (15 m/min) and recovery show that total lipid oxidation of a mouse increases during the first minute of exercise, but decreases to resting levels during steady-state exercise. Researchers could also selectively enrich the proteins in the body by feeding rodent diet infused with ^{13}C -L-Leucine, thereby allowing quantitative assessment of endogenous protein oxidation in real-time.



SPECIFICATIONS

	$\delta^{13}\text{C}$: 3.0‰ / 1.0‰ / 0.30‰
PRECISION	$\delta^{18}\text{O}$: 16‰ / 5‰ / 2‰
(1 SEC/10 SEC/100 SEC)	$^{12}\text{CO}_2$: 0.8 ppm / 0.30 ppm / 0.10 ppm
	H_2O : 200 ppm (60 sec) / 100 ppm (300 sec)
TOTAL UNCERTAINTY	< 1%
MEASUREMENT RANGE	CO_2 : 380 – 25000 ppm
(MEETS ALL SPECS)	H_2O : 4000 – 60000 ppm
OPERATIONAL RANGE	H_2O : 0 – 70000 ppm (non-condensing)
	CO_2 : 0 – 50000 ppm
MEASUREMENT RATES	User-selectable rates up to 1 Hz
	Ambient Humidity: non-condensing (0 – 100% RH)
SAMPLING CONDITIONS	Sample Temperature: -20 – 50 °C
	Operating Temperature: 0 – 45 °C
FITTINGS	Outlet (internal pump): 1/4 in.
	Inlet: 3/8 in.
OUTPUTS	Digital (RS-232), Ethernet, USB
POWER REQUIREMENTS	115/230 VAC, 50/60 Hz 66 W
DIMENSIONS (H x W x D)	18 cm (7") x 47 cm (18.5") x 36 cm (14")
WEIGHT	17 kg (37.5 lbs)

For more information on stable isotope labeling utilization, see:

McCue, M. D. (2011). "Tracking the oxidative and non-oxidative fates of isotopically labeled nutrients in animals." *BioScience* 61(3): 217-230.

Welch Jr, K. C., et al. (2016). "Combining respirometry with stable isotopes to investigate fuel use in animals." *Annals of the New York Academy of Sciences* 1365(1): 15-32.

McCue, M. D. and K. C. Welch Jr (2016). " ^{13}C -Breath testing in animals: Theory, applications, and future directions." *J Comp Physiol* 186B(3): 265-285.

ABOUT US

Sable Systems International designs and manufactures leading-edge gas, metabolic and behavioral measurement systems for calorimetry, respirometry, metabolic/behavioral phenotyping, and gas analysis. Our products enable the highest precision and resolution, optimum workflow and reliable performance – giving you utmost confidence in your results. By scientists, for scientists, Sable enables results that impact research and industry breakthroughs.



Sable Systems International
3840 N. Commerce Street
North Las Vegas, NV 89032, USA
TELEPHONE:
US: +1 800 330 0465 / + 1 702 269 4445
EMAIL: sales@sablesys.com

www.sablesys.com



Sable Systems Europe GmbH
Ostendstr. 25
D-12459 Berlin, Germany
TELEPHONE: +49 30 5304 1002
MOBILE: +49 176 2078 7008
FAX: +49 30 5304 1003